WHAT'S NEW IN VECTRA POLARIS 1.0?



Vectra Polaris scans both brightfield and fluorescent slides, and can take multispectral imagery of user-specified areas of interest, user-chosen TMA cores, or areas of interest determined through the use of inForm Tissue Finder.

Below, we've highlighted the new features of Vectra Polaris 1.0

- Slide Carriers
- Dashboard
- Protocol Editor
 - Opal Kit-Based Protocols (new in 1.0.8!)
 - Multispectral Whole Slide Scanning (new in 1.0.8!)
- Snapshots and Exposures
- Scanning Slides







Vectra Polaris uses slide carriers for touchless slide processing; up to 4 slides can sit in a single carrier.

The Vectra Polaris system can hold 20 carriers at a time.




HOME SCREEN

The Vectra Polaris home screen is organized by task.

The Dashboard provides the tools to confirm that your system is configured properly before scanning.

The Protocol Editor is a one-stop place for protocol creation and editing.

The slide scanning editor lets you easily set up your scans and walk away while Polaris scans.





The Vectra Polaris Dashboard provides tools to let you:

- See the available disk space.
- View and take brightfield and fluorescence references.

/iew Fluorescence Refere	ences	Disk Space Location: C:\Data\Jenga
Select Objective: 20x Scan Select Filter: Texas Red	Reference Information Available: Yes Last measured: Apr 20, 2016 13:11:15 Exposure Time: 251.17 (ms) Peak Counts: 181.65 Responsivity: 2.84	132 GB free Brightfield References View Acquire
	Export For Diagnostics	Ruorescence References View Acquire

Vectra Polaris




DASHBOARD

PROTOCOL CREATION

Users who upgrade their filters with the 1.0.8 software release are able to easily create protocols for Opal stained slides.

- You can scan whole slides multispectrally with Opal Polaris 7 and Opal 4 staining.
- A conventional whole slide scan with multispectral fields is available for the Opal 7 kit, as well as Opal 8 and 9 color.

Select the Opal option in the Staining section to create a protocol with the necessary filters, or you can select custom filters for imaging as well.





Protocol Editing in Vectra Polaris is a simple process:

- Select your imaging resolutions for the slide scan and multispectral regions.
- Set your filters, bands, and exposures in fluorescence.

PROTOCOL EDITING

👖 Vectra Polaris					-	×
Back Edit	t Protocol					-
New Loa	ad Save Protocol: T	estFL				
Overview Scan Rules						
Overview Scan Filter:	DAPI ~					
Multispectral Slide Scar	Settings	Multispectral Field Settin	76			
Pixel Resolution:	0.50µm (20x) ~		0.50µm (20x)	~		
Filter: DAPI / Opal570) / 690	Filter: DAPI / Opal 780				
DAPI	12.00 ms	DAPI MSI	DAPI	150.00 ms		
Opal 570	12.00 ms	Opal 780 MSI	Opal 780 or Cy7	150.00 ms		
Opal 690	12.00 ms	Filter: Opal 480 / Cy5				
Filter: Opal 480 / 620	/ 780	Opal 480 MSI	Opal 480	150.00 ms		
Opal 480	12.00 ms	Cy5 MSI	Opal 690 or Cy5	150.00 ms		
Opal 620	12.00 ms					
Opal 780	12.00 ms	Filter: FITC		150.00		
Filter: Opal 520		FITC	FITC or AF488	150.00 ms		
Opal 520	12.00 ms	Filter: CY3				
		Cy3	Cy3 or TRITC	150.00 ms		
Filter: Sample AF		Filter: Texas Red				
Sample AF	100.00 ms	Texas Red	Texas Red or	150.00 ms		
Scan Exposures	Select Scan Bands	Field Exposures	Select Field B	ands		
Advanced Settings						



PROTOCOL EDITING



Use the Advanced Settings button to access the Autofocus, TMA Sample finding, and Barcode settings in the protocol.

• If your samples are TMAs, check Sample is a TMA to use a tissue finding algorithm that is customized for the TMA format.

Sample Finding Options	
Sample is a TMA	

• If your slide labels have barcodes, you can store the decoded label contents within the .qptiff. Select the type of barcodes on your slide.

Label Barcode Reading	
No barcodes \sim	
No barcodes	
All 1D barcodes (autodet	ect type)
All 2D barcodes (autodet	ect type)
All barcode formats (auto	detect type)
Datamatrix only	
PDF 417 only	
QR Code only	
Aztec only	
Maxicode only	


170 µm (#1 or #1.5)		
	\sim	
Autofocus		
DAPI	\sim	For Slide Scans
DAPI MSI	\sim	For Fields
Saturation Protection		
Use Saturation Pr	rotection Fo	or Fields
Sample Finding Optio	ons	
Sample is a TMA		
Label Barcode Read	ing	
No barcodes	\sim	


SNAPSHOTS AND EXPOSURES



When setting exposures in fluorescence (or taking brightfield snapshots) you can take an overview to help locate relevant areas of tissue.

- Click the Take Overview button when a carrier is on stage to acquire an overview.
- A color overview is available for brightfield snapshots.
- An enhanced grayscale overview is available for fluorescence exposures.





Configure your scanning tasks for one or more carriers. You can apply the same settings to all slides, or different settings for individual slides.

Click Scan to start scanning!

You can pause a scan to unload completed slides and reload the system with new carriers to process.

_ Vectra Polaris \times Back Configure Enter Slide Scan Tasks IDs Θ Slot 1 🛸 slide2 Θ Slot 2 Θ \bigcirc \bigcirc Slot 3 Slot 4 Slot 5 e Θ Θ Θ







After configuring you scanning tasks, you can save the setup information for later use. Click Save Setup to save your selected studies, tasks, protocols, and slide IDs.

Click Load Setup to import previously saved work into the Scan Slides Editor.

You may need to place your carriers back in the hotel. See "Setting Up Scan Rules" in the manual for more information.



		Name	Notes	Studies	# Carriers	Permanent
6/15/2018 5:28:3	88 PM	PD-L1	520 titrations	180614 PDL1 520 titration and 180614 CD8 480 titrations	1	No
6/15/2018 5:27:1	2 PM	PDL1, CD8 titrations	CC, Opal titrations, tonsil	180614 PDL1 520 titration and 180614 CD8 480 titrations	2	No