

We've rebranded some of our products. *Phenoptics™* is now *Phenolmager™* and *Vectra® Polaris™* is now *Phenolmager™ HT*.

Opal Storage and Stability

■ LYOPHILIZED

Lyophilized Opal fluorophores (Opals) are extremely stable long-term when stored in the dark at -20°C. Performance should not be impacted if lyophilized Opal fluorophores have been exposed to room temperature (RT), if exposure was less than seven days.

■ RECONSTITUTED

Reconstitute all Opals except Opal 780 with 75ul Dimethylsulfoxide (DMSO). Reconstitute Opal Polaris 780 with 300uL double distilled water (ddH₂O). Reconstituted Opals should be stored at 4°C in the dark and must be used within three months when stored as directed.

■ DILUTED TO WORKING SOLUTION

Reconstituted Opals, with the exception of Opal Polaris 780, should be diluted in 1X Plus Amplification Diluent if staining manually or 1X Plus Automation Amplification Diluent if staining on automated platform. Opal Polaris 780 which is diluted with Opal Antibody Diluent/Block. Opal working solution must be stored at 4°C. Working solutions are good for 24-48 hours when stored properly. Prepare working solution for the volume needed, including automatic staining dead volume. Minimize exposure to light. DO NOT FREEZE working solution. Discard unused working solution.

OPAL RECONSTITUTION

■ BEFORE RECONSTITUTION

- Bring lyophilized Opal fluorophore and DMSO or ddH₂O to RT and ensure that DMSO is completely liquid before reconstitution.
 - DMSO freezes at 4°C. Ensure DMSO is free of crystals.
- Centrifuge lyophilizate for 10 seconds at 10 x g at RT.
- Dissolve lyophilizate in DMSO by gently vortexing.

■ AFTER RECONSTITUTION

- Centrifuge stock solution for 10 seconds at 10 x g at RT.
- Let Opal fluorophore stock stand 10-15 minutes at RT to ensure lyophilizate has completely dissolved.
- Vortex and centrifuge Opal fluorophore stock again for 10 seconds at 10 x g at RT.
- Prepare working aliquots of 5 - 10 ul each depending on desired use and wrap aliquot caps with parafilm.
 - DMSO is extremely hygroscopic and water could cause the dye to hydrolyze.
- Store in the dark at 4°C.
 - Use non-transparent caps OR wrap aliquots in aluminum foil.
- For subsequent preparation of working solution please ensure reconstituted Opals are completely dissolved before pipetting.
 - Hold in your hand and flick the tube a few times before gently vortexing.
 - Important: Centrifuge solution for 10 s at 10 x g at RT.
 - Open the cap and pipette the minimal amount required for a working solution.
- Minimize the number of freeze/thaw cycles by aliquoting reconstituted Opal fluorophore stock solution into microcentrifuge vials.

To learn more visit [AKOYABIO.COM](https://www.akoynbio.com) or email us at INFO@AKOYABIO.COM

For Research Use Only. Not for use in diagnostic procedures.

© 2022 Akoya Biosciences, Inc. All rights reserved. All trademarks are the property of Akoya Biosciences unless otherwise specified.

DN-00119 Rev A