

# Profiling spatial signatures in the microenvironment of non-small cell lung carcinomas in immunotherapy responsive and resistant disease cohorts




Ning Ma<sup>1</sup>, Aditya Pratapa<sup>1</sup>, James Monkman<sup>2</sup>, Najiba Mammadova<sup>1</sup>, Gabrielle Belz<sup>2</sup>, Ken O'Byrne<sup>3</sup>, Oliver Braubach<sup>1</sup>, Arutha Kulasinghe<sup>2</sup>




<sup>1</sup>Akoya Biosciences, Menlo Park, CA, USA; <sup>2</sup>The University of Queensland Diamantina Institute, The University of Queensland, Brisbane, QLD, Australia; <sup>3</sup>Princess Alexandra Hospital, Brisbane, QLD, Australia

## 1. Introduction


Immunotherapies have revolutionized cancer treatment. The immune contexture of the tumor microenvironment (TME) is an important factor in dictating how well a tumor responds to immune checkpoint inhibitors, but biomarker information alone is insufficient to design appropriate treatments. Instead, spatial classification of the immune contexture within the TME is foundational to addressing how immunological composition and status dictate patient responses. We used a combination of high-fidelity and high throughput spatial biomarker imaging to thoroughly characterize the TME of non-small cell lung carcinoma (NSCLC) tissue samples against outcomes for the immunotherapy arm of the Nivolumab trial. Tissues were sourced from cohorts of complete responders, partial responders, stable disease (not shown), and progressive disease patients. We first analyzed tissues using an ultra-high plex antibody panel of 46 antibodies with the PhenoCycler™-Fusion system, which provides whole-slide, single-cell resolution imaging data from FFPE samples. We then analyzed serial sections via whole-slide multiplex-immunohistochemistry (mIHC) using the Phenomager™ HT spatial imaging platform. Here, we used six-plex antibody panels – derived from the original 46-plex panel - aimed at T-cell profiling, immune exhaustion, and immuno-modulatory targets. Both analyses revealed similar spatial signatures across cohorts, which were all derived from the immunotherapy arm of the study.



## 2.1 PhenoCycler™-Fusion: Ultrahigh-plex Spatial Phenotyping




## 2.2 Phenomager™ HT: High throughput Spatial Phenotyping



## 3. Ultrahigh-plex Spatial Phenotyping of Non-small Cell Lung Cancer Tissue Cohorts

### A) 46-plex Whole-slide Imaging




46-plex Deep Immune Profiling Panel			
Immune Cell Core	Structural Module	Activity Module	Advanced Module
CD4	CD14	E-cadherin	PD-1*
CD68*	CD56	CD31	PD-L1*
CD20	CD19	Podoplanin	ICOS
CD11c		SMA	FOXP3*
CD8*		Vimentin	TIM3
HLA-DR		Collagen IV	LAG3
Ki67		CD34	IDO1
CD45RO	CD163	B-catenin	CD40
PanCK*	CD11b		CD40
CD3e	CD11b		HLA-E
CD44	MPO		IFNG
CD45	iNOS		VISTA
			CD15
			HLA-A


\*Same markers are in MOTIF™ Panel below

**A) PhenoCycler™-Fusion whole-slide 46-plex imaging** of non-small cell lung carcinoma tissue from different immunotherapy cohorts (biomarkers as indicated, scale bar, 2mm). The antibody panel on the left was designed for deep-phenotyping of immune cells. **B) Single-cell Spatial Phenotyping** of 539,158 cells in an NSCLC PR tissue. Cells are phenotyped via unsupervised clustering of biomarkers. The abundance of each cell type is quantified in the accompanying pie chart. **C) Cellular neighborhoods in NSCLC tissue.** Seven different neighborhoods (colored areas) contain different proportions of immune, cancer and epithelial cells. The heatmap shows the log fold change from each cell neighborhood compared to the average of the whole tissue, which indicates the respective frequencies of different cell types within each cellular neighborhood.

### B) Single Cell Spatial Phenotyping




### C) Cellular Neighborhoods




## 4. Multispectral Spatial Analysis Predicts Clinical Response to Immunotherapy

### A) High-throughput whole-slide immunofluorescence imaging




Serial sections from the same tissues as in Box 3 above were imaged with MOTIF™ on the Phenomager™ HT platform.

### B) Spatial analysis of NSCLC tissue cohorts



## 5. Correlation of Spatial Phenotyping Data with DSP Readouts



## 6. Conclusions and Outlook

- This study reports interim analysis for a multi-omic dataset from a single-agent Nivolumab arm from a clinical trial and the data are reliable across different methods, tissue cohorts and collection sites.
- Imaging data were quantified using different approaches and resulting data were consistent with prior findings. Our data demonstrated that high-plex analysis of the TME is likely to provide new insights into therapy resistance, and sensitivity.
- Identifying the benefit of treating patients with immune checkpoint inhibitor therapy prior to treatment will ensure better outcomes. Considering that tumor tissues are highly variable and change dynamically, there is an urgent need for the development of additional and more predictive biomarkers of response to IO therapy. Spatial phenotyping represents a solution for this unmet clinical need.