


We've rebranded some of our products. **CODEX®** is now **PhenoCycler™**.

Conjugating PhenoCycler Tags on Antibodies of Choice

1. Prepare Conjugation Reagents and Antibodies

- Purify antibodies if needed (see User Manual)
- Measure and calculate concentration of Antibodies using pre-set IgG settings of Nanodrop
- Calculate the volume of solution corresponding to **50 µg of antibody**.
- Retrieve the following reagents now:
 - Reduction Solutions 1 & 2
 - Filter Blocking Solution
- Retrieve the following reagents after 30-min incubation for step 4.1.5:
 - Conjugation Solution
 - Barcodes
- Retrieve the following reagents after 2-hour incubation for step 4.1.8:
 - Purification Solution
 - Antibody Storage Solution



2. Reduce Purified Antibody

- Label a 50kDa MWCO filter for each antibody.
- Add **500 µL of Filter Blocking Solution** to the top of each **50kDa MWCO filter**.
- Spin down at **12,000g** for **2 mins**.
- Remove all liquid; discard both the liquid left on the top of the column and the flow-through solution. Use a pipette if necessary to remove all liquid on top of the column. Add **50 µg** of the **antibody** in a volume **100 µL** or greater.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.
- Create the **Antibody Reduction Master Mix** based on the number of PhenoCycler antibody conjugates.

ANTIBODY REDUCTION MASTER MIX								
Number of Antibodies	1	2	3	4	5	6	7	8
Reduction Solution 1 [µL]	6.6	13.2	19.8	26.4	33	39.6	46.2	52.8
Reduction Solution 2 [µL]	275	550	825	1100	1375	1650	1925	2200
Total [µL]	281.6	563.2	844.8	1126.4	1378	1689.6	1971.2	2252.8

- Add **260 µL** of the **Antibody Reduction Master Mix** to the top of each filter unit. Vortex for 2-3 seconds.
- Incubate at RT for **30 mins**. Exceeding 30 mins will result in irreparable damage to antibodies.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.
- Add **450 µL** of **Conjugation Solution** to the top of each column.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.

Antibody stocks

Purified Antibodies in PBS or in a similar buffer. Antibodies should be free of carrier proteins, such as BSA, and other chemicals such as gelatin or glycerol.

Akoya Materials

- PhenoCycler Conjugation Kit
 - Reduction Solution 1
 - Reduction Solution 2
 - Filter Blocking Solution
 - Conjugation Solution
 - Purification Solution
 - Antibody Storage Solution
 - PhenoCycler Barcodes

Materials NOT Included in Kit

- Biologics:
 - Purified antibody(s)
- Consumables:
 - 50kDa MWCO filter
 - 1.5 mL screw-top sterile tube(s)
 - Nuclease-Free Molecular Biology Grade Water
 - PBS
 - 0.2 mL PCR tubes (for quality check, Section 4.2)
 - Bucket of ice for antibodies
- Instrumentation:
 - Centrifuge for 1.5 mL tubes
 - NanoDrop™ spectrophotometer
 - Vortex (Optional)

3. Prepare PhenoCycler Barcodes and Conjugate Antibody


- Resuspend each **Barcode** in **10 µL** of Molecular Biology Grade **Water** by pipetting up and down.
- Add **210 µL** of **Conjugation Solution** to each **suspended Barcode**. Mix by gentle pipetting. Set aside.
- Add the **PhenoCycler Barcode Solution** to the top of each filter. Close the lid and vortex for 3 seconds.
- Incubate for **2 hours** at **RT**.

4. Purify and Collect Tagged Antibodies

- Set aside **2-5 µL** of the purified solution for QC and troubleshooting.
- Spin down tubes at **12,000 g** for **8 mins**. Discard flow-through.
- Add **450 µL** of **Purification Solution** to the top of each column.
- Spin down tubes at 12,000 g for 8 mins. Discard the flow-through.
- Repeat steps c and d** two more times for **a total of three** purifications.
- After the **third centrifugation**, **discard** the flow-through solution.
- The top of the column will contain the remaining purified solution. The filter will contain the conjugated antibody solution.
- For each antibody, **label a new tube** and lid that can hold filter units.
- Add 100 µL of Antibody Storage Solution to each filter unit column.
- After it is labeled, place the new empty tube upside-down on top of the filter unit column.
- Invert the filter** for collection into the new collection tube.
- Spin solution down at **3,000g** for **2 mins**. The **final volume** in the tube should be **around 120 µL**.
- Transfer the solution to a sterile, screw-top tube for storage at 4°C for up to 1 year.
- Do not use these antibodies for tissue staining for at least 2 days; if used for staining sooner, you may experience high levels of background nuclear staining.

5. Validating Custom-Conjugated Antibodies

- Dilute** conjugated antibodies and the unconjugated antibody (control) to a **final volume** of **13 µL** in Nuclease free water.
- Add **5 µL** of NuPAGE™ LDS Sample Buffer to each sample.
- Add **2 µL** case NuPAGE™ Sample Reducing Agent to each sample.
- Denature at **95°C** in a dry bath for **10 mins**.
- Prepare buffer by adding **40 mL** of NuPAGE MOPS SDS Running Buffer and 760 mL of NuPAGE MOPS SDS Running Buffer.
- Prepare gel according to manufactures instructions.
- Pour buffer into gel tank.
- Add **5 µL** of a pre-stained protein standard 3.5-260 kDa to the gel.
- Add **20 µL** of antibody to one well each.
- Run the gel **at 200 V** for **30-40 mins** until completion.
- Place in a container filled with ddH₂O.
- Microwave** the gel until the **first bubbles form**.
- Stain the gel with Novex SimplyBlue™ SafeStain.
- Microwave** the gel again until the first **bubbles form**.
- Place the gel in a shaker for **10 mins**.
- Wash the gel with ddH₂O and leave it on the shaker until bands are visible.



From left to right: the first column shows the protein standard, columns from the second to the fifth show bands of successfully barcode-conjugated antibodies.

The last column shows the heavy and light chain bands from a control, unconjugated antibody.

Antibody

- 5 µL of each conjugated antibody from section
- 4.1.9a.
- 1 µg (usually corresponding to 2 µL) of unconjugated antibody to be used as control.

Materials NOT Included in Kit

- Reagents:
 - NuPAGE™ LDS Sample Buffer (4X)
 - NuPAGE™ Sample Reducing Agent (10X)
 - NuPAGE™ 4-12% Bis-Tris Protein Gels
 - Novex™ Sharp Pre-Stained Protein Standard -3.5260 kDa
 - XCell SureLock™ Mini-Cell Electrophoresis System
 - NuPAGE™ MOPS SDS Running Buffer (20X)
 - Novex™ SimplyBlue™ SafeStain
 - ddH₂O
 - Nuclease-free water
- Instrumentation:
 - 95°C dry bath
 - Microwave