


6768: The Potential Predictive Role of Spatial Phenotyping in Non-small Cell Lung Cancer (NSCLC)

Ning Ma¹, Aditya Pratapa¹, James Monkman², Ken O'Byrne³, Oliver Braubach¹, Arutha Kulasinghe²

¹Akoya Biosciences, Inc, Marlborough, Massachusetts, USA; ²University of Queensland, Brisbane, QLD, AUS; ³Princess Alexandra Hospital, Brisbane, QLD, AUS

1. Introduction


Lung cancers are the leading cause of cancer-related deaths with a 5-year survival of only ~20%. Whilst immunotherapies have led to durable and prolonged survival, only a subset of patients remains responsive. Additional biomarkers are thus needed to better predict if patients will respond or develop resistance against immune checkpoint inhibitor (ICI) therapies. Spatial phenotyping of the tumor microenvironment (TME) is now recognized as a predictive proxy for ICI therapy outcomes. We phenotyped pretreatment biopsies from non-small cell lung cancer (NSCLC) patients enrolled in a single-agent Nivolumab clinical trial. We first performed 57-plex whole-slide Single Cell Spatial Phenotyping on the PhenoCycler®-Fusion platform. Our analyses revealed high phenotypic diversity in the TME of responding and non-responding patients. We then deployed customizable **PhenoCode™ Signature Panels (PSP)** for high-throughput immune profile (IP: CD3ε/CD8/CD20/CD68/PanCK + CD4 add-in) and immuno-contexture (IC: CD8/CD68/PD-L1/FoxP3/PanCK + PD-1 add-in) analyses on the PhenImager® HT. The PSP panels were deployed on 37 biopsies and afforded comprehensive evaluation of patient cohorts. Our PSP data revealed no difference in the immune cell makeup in responder vs. non-responder TMEs. However, we discovered multiple quantifiable and statistically significant **Spatial Signatures** that appear to be predictive of a positive treatment outcome.



2. Enabling IO Discovery to Translation with Akoya's Unified Solutions




PhenoCode Signature Panels use Akoya's barcode-based antibody labeling chemistry and are validated for the PhenImager® HT workflow. Featuring a flexible design component that allows for the easy integration of a novel checkpoint or immune cell marker, these panels offer 3-fold faster assay development and optimization times when compared to other custom 6-plex panels.




PhenoCode Signature Panels use Akoya's barcode-based antibody labeling chemistry and are validated for the PhenImager® HT workflow. Featuring a flexible design component that allows for the easy integration of a novel checkpoint or immune cell marker, these panels offer 3-fold faster assay development and optimization times when compared to other custom 6-plex panels.

3. Ultrahigh-plex & High-throughput Single-Cell Spatial Phenotyping Analyses of NSCLC Cohorts


57-plex Spatial Phenotyping Identified 11 Immune Cell Populations in Responder vs. Non-responder



6-plex Spatial Phenotyping Validates Diverse NSCLC Immune Landscape but Rules out Differences in Immune Composition




Immuno-Contexture PSP, Spatial Phenotyping with PSP-IC, and Quantification of Immune Activity




4. Cell Neighborhood Analyses Indicate Distinct Spatial Biology of Responder vs. Non-Responder TME

NSCLC Biopsies Contain 7 Distinct Cellular Neighborhoods



5. A Predictive Spatial Signature for NSCLC Treatment Outcomes



The **SpatialScore** is the ratio of the physical distance between CD8+ T-cells and the nearest tumor cell, relative to its nearest Treg (top) or Macrophage (bottom). Both, Treg and Macrophages can exert immunosuppressive effects on CD8+ T-cells, which then results in a higher SpatialScore. Indeed, the pooled SpatialScore from pretreatment biopsies were significantly higher in the Non-responder (NR) condition when compared to the responder (R) condition. A high spatial score can thus be interpreted as: higher CD8+ T-cell suppression - lower anti-tumor activity - lower survival rate (see Kaplan Meier curve in Introduction).

6. Conclusions and Outlook

- This study amounts to a uniquely comprehensive Single Cell Spatial Phenotyping analysis of **pretreatment** NSCLC biopsies from a single-agent Nivolumab clinical trial.
- Our data illustrate the diverse immune microenvironment of NSCLC and indicate that immune cell quantification is insufficient to stratify patient cohorts. The identification of unique spatial signature such as a spatial score may be developed to predict patient response.
- Single-cell Spatial Phenotyping with Akoya's unified PhenoCycler-Fusion and PhenImager solution along with PhenoCode Discovery and Signature panels enables the discovery of multiple quantifiable and statistically significant Spatial Signatures in Responder vs. Non-responder cohorts.
- This study shows how Akoya's solutions are uniquely positioned to enable discovery to translational workflows thereby accelerating the development of clinically relevant and highly predictive spatial signatures.**

¹Phillips et al. (2021). Immune cell topology predicts response to PD-1 blockade in cutaneous T cell lymphoma. *Nature Comm.*