


# 6872: A Spatio-Temporal Approach to Mapping the Dynamics of Cutaneous Squamous Cell Carcinoma Progression and Immunotherapy Response: A Journey through TiME

Niyati Jhaveri<sup>1</sup>, Bassem B. Cheikh<sup>1</sup>, Dmytro Klymyshyn<sup>1</sup>, Ning Ma<sup>1</sup>, Aditya Pratapa<sup>1</sup>, James DeRosa<sup>1</sup>, Ritu Mihani<sup>1</sup>, James Monkman<sup>2</sup>, Nadine Nelson<sup>3</sup>, Arutha Kulasinghe<sup>2</sup>  
<sup>1</sup>Akoya Biosciences, Marlborough, California, USA; <sup>2</sup>University of Queensland, Brisbane, AU; <sup>3</sup>Abcam, Cambridge, UK

## 1. Cutaneous Squamous Cell Carcinoma


**Cutaneous Squamous Cell Carcinoma (cSCC)**, arising from proliferation of malignant epidermal keratinocytes, is the *second most common non-melanoma skin cancer*<sup>1</sup>. Tumor development is a gradual process characterized by a high *mutational burden* and an *immunosuppressive* tumor immune microenvironment (TiME). Though **immunotherapy** is a promising solution, and prognoses are favorable in most cases with a 5-year survival  $\geq 90\%$ , cSCC accounts for 75% of all skin-cancer related deaths excluding melanoma<sup>1</sup>.



<sup>1</sup>Fania et al., Biomedicines, 2021 Feb 9;9(2):171. Illustration: Centers for Disease Control and Prevention (CDC)

## 2. PhenoCycler-Fusion 2.0 Technology

The **PhenoCycler®-Fusion 2.0 (PCF)** technology is the fastest whole-slide spatial biology platform that enables simultaneous detection of 100+ biomarkers by combining automated fluidics and iterative imaging of DNA-tagged antibodies. Leveraging the high plex, high throughput single-cell resolution of PCF, we sought to identify the **spatio-temporal changes** in the cSCC microenvironment. In this interim analysis, we profiled two patients before immunotherapy with Cemiplimib and 6-12 months following treatment.




### PhenoCode™ Discovery Panels


<b>Cell Profiling Core</b>	CD4, CD68, CD20, CD11c, CD8, HLA-DR, Ki67, CD45RO, CD3e, PanCK, CD44, CD45, HLA-A, CD14, CD56 + CD163, CD11b, CD57
<b>Lymphocyte Profiling Module</b>	CD21, FOXP3, Granzyme B, CD79a, TCF-1, CD38, TOX, T-bet, CD107a, CD39
<b>Immune Activation and Proliferation Module</b>	PD-1, PD-L1, IDO1, LAG3, VISTA, IFNG, HLA-E, CD40, ICOS, PCNA
<b>Tissue Architecture Module</b>	E-cadherin, CD31, Podoplanin, SMA, Vimentin, Collagen IV, CD34, $\beta$ -catenin1, $\beta$ -actin, Caveolin
<b>Custom Metabolic Module</b>	ASCT2, CPT1A, HK1, LDHA, IDH2, GLUT1, pNRF2, ATP5A, Citrate Synthase, HIF1a, G6PD, SDHA

## 3. Spatio-Temporal Analysis Reveals Differential Tumor-Immune Phenotypes


### 3.1 Whole-Slide Single-Cell Spatial Phenotyping Identifies 12 Distinct Cellular Phenotypes in cSCC Tissues



### 3.2 Spatio-Temporal Analysis Reveals a “Responder” Phenotype with Tumor Regression and M1 Polarization




### 3.3 Spatio-Temporal Analysis Reveals a “Non-Responder” Phenotype with Tumor Burden and M2 Polarization




## 4. Spatial Analyses Reveals Regional Heterogeneity

### 4.1 Spatial Proximity Analyses Reveals Temporal Differences in Cellular Distribution



### 4.2 Tissue Compartmental Analysis Highlights Regional Heterogeneity



## 5. The Power of Spatio-temporal Mapping at Scale

**Cancer** development is a dynamic process, characterized by cellular and molecular changes in the tissue microenvironment that contribute to sustained proliferation, immune evasion, and resistance. Since cSCC biopsies can be collected non-invasively, longitudinal studies tracking spatio-temporal changes can follow the course of the tumor in the same patient. The PCF 2.0 platform is ideally suited for such studies with its high plex, high throughput single-cell resolution, seamless scaling, and unparalleled speed.



Download Poster