


## 1. Introduction


The transition from ultrahigh-plex discovery experiments in multiplexed imaging to targeted, high-throughput translation and thereby clinical research studies presents a significant challenge. Successful translation hinges on the harmonization of staining, imaging, and data analysis technologies across platforms. This work aims to present a streamlined workflow for the transition, leveraging the capabilities of the PhenoCycler®-Fusion 2.0 and Phenolmager® HT platforms.

We demonstrate a two-step approach utilizing PhenoCode™ panels: 1. **Deep Protein Spatial Phenotyping:** For this step, tumor microarrays (TMAs) comprised of human formalin-fixed, paraffin-embedded (FFPE) cancer tissues were phenotyped using the ultra-high-plex PhenoCode Discovery panels on the PhenoCycler-Fusion 2.0 platform. This step allows for extensive characterization of the tissue microenvironment at the protein level. 2. **Targeted High-Throughput Analysis:** Following discovery, PhenoCode Signature panels were employed on the Phenolmager HT 2.0 platform. This platform enables high-throughput analysis of protein targets identified during the discovery phase, facilitating the translation of these findings into targeted translational and clinical research studies.

## 2. Discovery to Clinical Research Workflow




## 3. Ultrahigh-Plex Spatial Phenotyping of Breast Cancer Tissues




**Single-Cell Spatial Phenotyping of Proteins in Breast Cancer Tissues using 32-plex Discovery Panel reveals distinct cell populations.** Three types of breast cancer TMA cores: triple negative breast cancer, invasive breast carcinoma, and invasive ductal carcinoma were imaged and analyzed. Heatmap (A) corresponding to average expression of cell-type defining markers revealed 13 distinct cellular populations identified by unsupervised clustering. (B) Representative images from invasive ductal carcinoma. The difference in the immune cell composition and status in the invasive ductal carcinoma were further analyzed using 6-plex PhenoCode Signature Panels which are ideally suited for automated high-throughput studies.

## 4. Quantitative Spatial Profiling of Breast Cancer Subtypes using 6-plex Multiplex Imaging and Phenotyping Analysis


### 4.1. Assessing immune cells infiltration in the tumor




### 4.2. Determining major immune cell types in the tumor



### 4.3. Analysis of T cell status within the tumor




## 4.4. Assessing macrophage polarization in the tumor






## 5. Assessment of Spatial Relationships in Breast Cancer Subtypes

### 5.1. Nearest Neighborhood Analysis




### 5.2. Spatial Score Analysis



**5.1. Nearest Neighbor (or spatial proximity) analysis could indicate prognosis in breast cancer.** The PSP Activated TIL panel + CD4 was used to determine the distance between activated cytotoxic T cells (CD8+/Granzyme B+) and proliferating tumor cells (PanCK+/Ki67+). Representative images for invasive breast and ductal carcinoma are shown (A & B). Shorter distance between cytotoxic T cells and tumor cells seen in invasive breast carcinoma (C) versus invasive ductal carcinoma (D) could indicate better prognosis as suggested in previous studies.

**5.2. Spatial score indicates T cell suppression in invasive ductal carcinoma.** The PSP Immuno-Contexture panel + CD4 was used to assess the distances of tumor cells and Tregs to CD4+ T Cells. Representative images of invasive breast (A) and ductal (B) carcinoma are shown. The per cell (C) spatial score (see below for spatial score details) were assessed showing a higher score in invasive ductal carcinoma samples.



## 6. Unlocking the Power of PhenoCode Panels for Ultrahigh-plex Discovery to Targeted High-throughput Studies

PhenoCode panels, powered by Akoya's novel barcode chemistry, provide an off-the-shelf solution for immunofluorescence staining that requires minimal optimization. Our study showcases how ultrahigh-plex profiling using PhenoCode Discovery panels can be followed by a more targeted analysis using PhenoCode Signature panels. Here, we have shown how the PhenoCode Signature panels can be leveraged for rapid immune profiling and immune status assessment of the TME and to perform spatial analysis to gain key biological insights that can be used for the development of spatial signature that can more reliably predict immunotherapy response. [akoyabio.com/phenocode/](http://akoyabio.com/phenocode/)

