

CODEX Community Antibody List

CODEX Community Antibodies: Antibodies that CODEX users have conjugated to CODEX Barcodes and used successfully in their experiments.

CODEX User: Simon Knott, Kenneth Gouin and Catherine Oh, Department of Biomedical Sciences, Cedars Sinai Medical Center, CA, USA			
Target	Clone	Vendor	Tissue
CD27	LG.3A10	Thermo Fisher	Mouse bone marrow (FF)
CD34	EP373Y	Abcam	Mouse bone marrow (FF)
CD47	Miap301	Thermo Fisher	Mouse bone marrow (FF)
CD62L	MEL-14	Thermo Fisher	Mouse bone marrow (FF)
CD117	2B8	Thermo Fisher	Mouse bone marrow (FF)
CD152/CTLA4	UC10-4B9	BioLegend	Mouse bone marrow (FF)
CD144/VE-Cadherin	BV13	Thermo Fisher	Mouse bone marrow (FF)
CD206	MR5D3	Thermo Fisher	Mouse bone marrow (FF)
TCR Gamma/Delta	GL-3	LSBio	Mouse bone marrow (FF)

CODEX User: Diane Saunders, Powers & Brissova Research Group, Vanderbilt University Medical Center, TN, USA			
Target	Clone	Vendor	Tissue
Arginase 1	Polyclonal (Cat # NBPI-32731)	Novus	Human pancreas (FF)
C-peptide	C-PEP-01	Thermo Fisher	Human pancreas (FF)
CD36	5-271	BioLegend	Human pancreas (FF)
CD86	BU63	Thermo Fisher	Human pancreas (FF)
CD163	GHI/61	BioLegend	Human pancreas (FF)
CD206 (MMR)	15-2	BioLegend	Human pancreas (FF)
Chromagranin A	LK2H10+PHE5+CGA/414	Novus	Human pancreas (FF)
Ghrelin	883622	R&D Systems	Human pancreas (FF)
Glucagon	K79bB10	Abcam	Human pancreas (FF)
Iba1	EPRI6589	Abcam	Human pancreas (FF)
iNOS	4E5	Thermo Fisher	Human pancreas (FF)

NTPDase3	B3S-C3	ectonucleotidases antibodies	Human pancreas (FF)
Pancreatic Polypeptide/PP	548416	R&D Systems	Human pancreas (FF)
Somatostatin	7G5	Novus	Human pancreas (FF)

CODEX User: Lauren Abell and Kimberly Smythe, Houghton Lab, Immunopathology Team, Fred Hutchinson Cancer Center, WA, USA

Target	Clone	Vendor	Tissue
CD10	SP179	Abcam	Human tonsil (FFPE)
CD11b	SP331	Abcam	Human tonsil (FFPE)
CD123	6H6	Invitrogen	Human tonsil (FFPE)
CD14	SP192	Abcam	Human tonsil (FFPE)
CD163	EPR19518	Abcam	Human tonsil (FFPE)
CLEC9A	EPR22324	Abcam	Human tonsil (FFPE)
CXCR2	EPR22301-103	Abcam	Human tonsil (FFPE)
DR5	EPR22276	Abcam	Human tonsil (FFPE)
FoxP3	236A/E7	Invitrogen	Human tonsil (FFPE)
HLA-DR	ab259258	Abcam	Human tonsil (FFPE)
LOX-1	EPR20750	Abcam	Human tonsil (FFPE)
NSE	EPR3377	Abcam	Human tonsil (FFPE)
PD-1	EPR4877(2)	Abcam	Human tonsil (FFPE)
PD-L1	CAL10	Abcam	Human tonsil (FFPE)