

## CODEX Antibody List: Human Fresh Frozen Targets

**CODEX Inventoried Antibodies Kit:** Fully validated antibodies conjugated to a distinct and activated CODEX Barcode, with the Reporter (Barcode-dye conjugate) complement for use in CODEX runs. Each kit contains one tube of Barcoded Antibody and one tube of the corresponding Reporter, sufficient volume for 25 CODEX runs (25 tests) at the standard staining conditions.

Akoya PN	CODEX Antibodies
4250005	CD2-BX002 (RPA-2.10)—Atto 550-RX002
4350008	CD3-BX015 (UCHT1)—Cy5-RX015
4350010	CD4-BX021 (SK3)—Cy5-RX021
4150004	CD8-BX004 (SK1)—Alexa Fluor™ 488-RX004
4450004	CD8-BX004 (SK1)—Alexa Fluor™ 750-RX004
4350012	CD11c-BX027 (S-HCL-3)—Cy5-RX027
4350003	CD19-BX003 (HIB19)—Cy5-RX003
4150009	CD21-BX013 (Bu32)—Alexa Fluor™ 488-RX013
4450009	CD21-BX013 (Bu32)—Alexa Fluor™ 750-RX013
4250009	CD31-BX032 (WM59)—Atto 550-RX032
4250020	CD34-BX035 (561)—Atto 550-RX035
4150007	CD38-BX007 (HB-7)—Alexa Fluor™ 488-RX007
4450007	CD38-BX007 (HB-7)—Alexa Fluor™ 750-RX007
4150003	CD45-BX001 (HI30)—Alexa Fluor™ 488-RX001
4450003	CD45-BX001 (HI30)—Alexa Fluor™ 750-RX001
4250023	CD45RO-BX017 (UCHL1)—Atto 550-RX017
4350007	CD49f-BX033 (GoH3)—Cy5-RX033
4250022	CD69-BX041 (FN50)—Atto 550-RX041
4150021	CD90-BX022 (5E10)—Alexa Fluor™ 488-RX022
4450021	CD90-BX022 (5E10)—Alexa Fluor™ 750-RX022
4250008	CD104-BX005 (58XB4)—Atto 550-RX005
4150008	CD138-BX010 (MI15)—Alexa Fluor™ 488-RX010
4450008	CD138-BX010 (MI15)—Alexa Fluor™ 750-RX010
4250013	CD278-BX017 (C398.4A)—Atto 550-RX017
4250010	CD279-BX014 (EH12.2H7)—Atto 550-RX014
4250006	HLA-DR-BX026 (L243)—Atto 550-RX026
4250019	Ki67-BX047 (B56)—Atto 550-RX047
4150020	Pan-Cytokeratin-BX019 (AE-1/AE-3)—Alexa Fluor™ 488-RX019
4450020	Pan-Cytokeratin-BX019 (AE-1/AE-3)—Alexa Fluor™ 750-RX019

4250004

Podoplanin-BX023 (NC-08)—Atto 550-RX023