

## CODEX Antibody List: Mouse Fresh Frozen Targets

**CODEX Inventoried Antibodies Kit:** Fully validated antibodies conjugated to a distinct and activated CODEX Barcode, with the Reporter (Barcode-dye conjugate) complement for use in CODEX runs. Each kit contains one tube of Barcoded Antibody and one tube of the corresponding Reporter, sufficient volume for 25 CODEX runs (25 tests) at the standard staining conditions.

Akoya PN	CODEX Antibodies
4350014	CD3-BX021 (17A2)—Cy5-RX021
4250016	CD4-BX026 (RM4-5)—Atto 550-RX026
4250007	CD5-BX017 (53-7.3)—Atto 550-RX017
4250017	CD8a-BX029 (53-6.7)—Atto 550-RX029
4150015	CD11b-BX025 (M1/70)—Alexa Fluor™ 488-RX025
4450015	CD11b-BX025 (M1/70)—Alexa Fluor™ 750-RX025
4350013	CD11c-BX030 (N418)—Cy5-RX030
4250014	CD19-BX020 (6D5)—Atto 550-RX020
4250015	CD21/35-BX023 (7E9)—Atto 550-RX023
4150014	CD24-BX022 (M1/69)—Alexa Fluor™ 488-RX022
4450014	CD24-BX022 (M1/69)—Alexa Fluor™ 750-RX022
4250001	CD31-BX002 (MEC13.3)—Atto 550-RX002
4150013	CD38-BX019 (90)—Alexa Fluor™ 488-RX019
4450013	CD38-BX019 (90)—Alexa Fluor™ 750-RX019
4250002	CD44-BX005 (IM7)—Atto 550-RX005
4150002	CD45-BX007 (30-F11)—Alexa Fluor™ 488-RX007
4450002	CD45-BX007 (30-F11)—Alexa Fluor™ 750-RX007
4150006	CD45R/B220-BX010 (Ra3-6B2)—Alexa Fluor™ 488-RX010
4450006	CD45R/B220-BX010 (Ra3-6B2)—Alexa Fluor™ 750-RX010
4350007	CD49f-BX033 (GoH3)—Cy5-RX033
4350016	CD71-BX027 (RI7217)—Cy5-RX027
4450001	CD90.2-BX001 (30-H12)—Alexa Fluor™ 750-RX001
4150001	CD90.2-BX001 (30-H12)—Alexa Fluor™ 488-RX001
4350005	CD169-BX015 (3D6.112)—Cy5-RX015
4150012	IgD-BX016 (11-26c.2a)—Alexa Fluor™ 488-RX016
4450012	IgD-BX016 (11-26c.2a)—Alexa Fluor™ 750-RX016
4150011	IgM-BX013 (RMM-1)—Alexa Fluor™ 488-RX013
4450011	IgM-BX013 (RMM-1)—Alexa Fluor™ 750-RX013
4250019	Ki67-BX047 (B56)—Atto 550-RX047

4350015	Ly6g-BX024 (IA8)—Cy5-RX024
4250003	MHC II-BX014 (M5/114.15.2)—Atto 550-RX014
4350006	TCRb-BX003 (H57-597 )—Cy5-RX003
4150005	Ter119-BX004 (TER-119)—Alexa Fluor™ 488-RX004
4450005	Ter119-BX004 (TER-119)—Alexa Fluor™ 750-RX004